

Laboratory Test Updates

(effective on September 14, 2026)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8614005	C-Reactive Protein (micro-CRP) (STAT)	X									
2	C8230803	Procalcitonin (STAT)	X									
3	C8388001	B-Natriuretic Peptide (NT-proBNP) (STAT)	X									
4	C8530201	Protein C, Antigen		X			X					
5	C8744901	Aspergillus Galactomannan Antigen				X						

1	C-Reactive Protein (micro-CRP) (STAT)		
	Test Code : 8614005 - Laboratory Services		
	Methodology : sandwich immunodetection - Performed : Mon-Sat 14:00 - Reported : Mon-Sat 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Capillary Blood	Capillary Tube	2-8°C

2	Procalcitonin (STAT)		
	Test Code : 8230803 - Laboratory Services		
	Methodology : FEIA - Performed : Mon-Sat 08:00 - Reported : Mon-Sat 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL
			2-8°C

3	B-Natriuretic Peptide (NT-proBNP) (STAT)		
	Test Code : 8388001 - Laboratory Services		
	Methodology : FLUOROMETRY - Performed : Mon-Sat 08:00 - Reported : Mon-Sat 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	2mL
			2-8°C

4	Protein C, Antigen			
	Test Code : 8530201 - Laboratory Services			
	Methodology : IMMUNOTURBIDIMETRIC - Performed : Tue 08:00 - Reported : Thu 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (Citrate), Frozen	Blue top (Sodium citrate)	0,5mL	<0°C
	Reference Ranges 			

Normal Values				
Gender	Age Range	Lower Value	Upper Value	Unit
	0-4 DAY	17.0	53.0	%
	4-29 DAY	20.0	64.0	%
	29-89 DAY	21.0	65.0	%
	89-179 DAY	28.0	80.0	%
	179 DAY - 1 YEAR	37.0	81.0	%
	1-5 YEAR	40.0	92.0	%
	5-10 YEAR	45.0	93.0	%
	10-120 YEAR	67.0	125.0	%

5

Aspergillus Galactomannan Antigen

Test Code : **8744901** - Laboratory Services

Methodology : **CLIA** - Performed : **Mon-Thu 08:00 Fri 08:00** - Reported : **Tue-Fri 18:00 Mon 18:00**

Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Bronchoalveolar Lavage	Biologic Fluid Tube	1mL	2-8°C
Serum	Separator Tube	1mL	2-8°C
Plasma (EDTA)	Lavender top (EDTA)	1mL	2-8°C