

Laboratory Test Updates

(effective on June 06, 2022)

10	C8392409	Fetal Blood (Cordocentesis)_ Prenatal Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis	X								
11	C8392100	Amniotic Fluid_TRIO Prenatal Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis	X								
12	C8392163	CVS_TRIO Prenatal Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis	X								
13	C8392163	Bartonella henselae Antibody, IgG		X	X						

Culture, Aerobic (All Sample Types)Test Code : **8707501** - Laboratory Services1 Method : **C/S** - Performed : **Mon-Sat 12:00** - Reported : **3 days 17:00**

Sample Name	Container	Amount	Transfer
Skin Material	Call the Laboratory		2-8°C
Surgical Material	Call the Laboratory		2-8°C
Pleural Fluid	Biologic Fluid Tube	1mL	2-8°C
Biologic Fluid	Biologic Fluid Tube	1mL	2-8°C
Nail Scraps	Plastic Collection Cup		2-8°C
Cervical Swab	Biologic Fluid Tube		2-8°C
Cord Blood	Plain Tube		2-8°C
Pericardial Fluid	Biologic Fluid Tube	1mL	2-8°C
Ear Discharge	Biologic Fluid Tube		2-8°C
Tissue Sample	Plastic Collection Cup		2-8°C
Rectal Swab	Swap		2-8°C
Nasopharyngeal Swab	Swap		2-8°C
Nasal Discharge	Swap		2-8°C
Urethral Discharge Material	Swap		2-8°C
Peritoneal Fluid	Biologic Fluid Tube	1mL	2-8°C
Cord	Plastic Collection Cup		2-8°C
Gastric Juice	Biologic Fluid Tube	1mL	2-8°C
Ophthalmic Discharge	Biologic Fluid Tube	0,5mL	2-8°C
Biopsy Material	Plastic Collection Cup		2-8°C
Joint Fluid	Biologic Fluid Tube	1mL	2-8°C
Other	Call the Laboratory		2-8°C
Cystic Fluid	Sterile Tube		2-8°C
Bone Marrow	Biologic Fluid Tube		2-8°C
Vulva	Swap		2-8°C
Gastric Fluid, Fasting	Biologic Fluid Tube	5mL	2-8°C
Prostatic fluid	Biologic Fluid Tube		2-8°C

Culture, Fungi

Test Code : 8710201 - Laboratory Services

2

Method : C/I/S - Performed : Mon-Sat 12:00 - Reported : 21 days 17:00

Sample Name	Container	Amount	Transfer
Nail Scraps	Plastic Collection Cup		18-24°C
Urine, Sterile	Plastic Collection Cup		2-8°C
Throat Swab	Swap		2-8°C
Skin Scraps	Plastic Collection Cup		18-24°C
Vaginal Discharge	Swap		2-8°C
Hair	Plastic Collection Cup		18-24°C
Other	Call the Laboratory		2-8°C
Whole Blood	Blood Culture Bottle	10mL	18-24°C
Sputum	Plastic Collection Cup		2-8°C
Peripheral Stem Cell	Sterile Glass Tube	1mL	18-24°C
Bronchoalveolar Lavage	Biologic Fluid Tube	1mL	2-8°C

Amniotic Fluid_ Prenatal Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis

Test Code : 8392275 - Genetic Diagnosis Tests

7

Method : HIGH RESOLUTION SNP ARRAY - Performed : Mon-Sat - Reported : 2 weeks

Sample Name	Container	Amount	Transfer
Amniotic fluid and whole blood	Injector and Lavender (EDTA)	20mL	18-24°C
Whole Blood	Lavender (EDTA)	3mL	18-24°C

CVS_Prenatal Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis

Test Code : 8392274 - Genetic Diagnosis Tests

8

Method : HIGH RESOLUTION SNP ARRAY - Performed : Mon-Sat - Reported : 2 weeks

Sample Name	Container	Amount	Transfer
CVS and Whole Blood	Transport Medyum ve EDTA(mor kapak)	3villus	18-24°C
Whole Blood	Lavender (EDTA)	3mL	18-24°C

	Abortion Material_Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis Test Code : 8392124 - Genetic Diagnosis Tests			
Method : HIGH RESOLUTION SNP ARRAY - Performed : Mon-Sat - Reported : 2 weeks				
9	Sample Name	Container	Amount	Transfer
	Abortion Material and whole blood	Transport Medyum ve EDTA(mor kapak)	1Piece	18-24°C

	Fetal Blood (Cordocentesis)_ Prenatal Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis Test Code : 8392409 - Genetic Diagnosis Tests			
Method : HIGH RESOLUTION SNP ARRAY - Performed : Mon-Sat - Reported : 2 weeks				
10	Sample Name	Container	Amount	Transfer
	Cord Blood	Lavender (EDTA) + Green (Heparin)	3mL	18-24°C
	Whole Blood	Lavender (EDTA)	3mL	18-24°C

	Amniotic Fluid_TRIO Prenatal Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis Test Code : 8392100 - Genetic Diagnosis Tests			
Method : HIGH RESOLUTION SNP ARRAY - Performed : Mon-Sat - Reported : 2 weeks				
11	Sample Name	Container	Amount	Transfer
	Amniotic fluid and whole blood	Injector and Lavender (EDTA)	20mL	18-24°C
	Whole Blood	Lavender (EDTA)	3mL	18-24°C

CVS_TRIO Prenatal Microarray + Chromosome Analysis + Rapid Aneuploidy Analysis
Test Code : **8392163** - Genetic Diagnosis Tests

Method : **HIGH RESOLUTION SNP ARRAY** - Performed : **Mon-Sat** - Reported : **2 weeks**

12

Sample Name	Container	Amount	Transfer
CVS and Whole Blood	Transport Medyum ve EDTA(mor kapak)	5villus	18-24°C
Whole Blood	Lavender (EDTA)	3mL	18-24°C

Bartonella henselae Antibody, IgG

Test Code : **8660908** - Laboratory Services

Method : **CLIA** - Performed : **Mon-Thu 08:00 Fri 08:00** - Reported : **Tue-Fri 18:00 Mon 18:00**

13

Sample Name	Container	Amount	Transfer
Serum	Separator Tube	1mL	2-8°C